


SURGICAL
TECHNIQUE GUIDE

INCLUDING



KING KONG
PLATING SYSTEM 2.5

Made in Germany

EN Working with AGOMED means working with a company committed to excellence. Our products are German engineered, biocompatible, and cutting-edge. Our team members collectively bring decades of medical technology experience to the table.

We're reliable. We're flexible. We're inventive.

AGOMED collaborates with distributors, hospitals, and doctors worldwide to create and develop exacting, state-of-the-art solutions for complex surgical problems. We work closely with specialists to insure that our implant systems for Traumatology/Orthopaedics of upper and lower extremities improve the patient's quality of life. Patient safety is always our number one priority.

DE Mit AGOMED zu arbeiten, bedeutet mit einem Unternehmen zu arbeiten, das sich zur Exzellenz und Hochleistung verpflichtet.

Unsere Produkte werden von deutschen Ingenieuren entwickelt, sind biokompatibel, innovativ und auf dem neuesten Stand der Technik. Unser Team hat jahrzehntelange Erfahrung im Bereich der Medizintechnik.

Wir sind zuverlässig. Wir sind flexibel. Wir sind innovativ.

AGOMED arbeitet weltweit mit Unternehmen, Krankenhäusern und Ärzten zusammen, um anspruchsvolle, hochmoderne Lösungen für komplexe chirurgische Probleme zu entwickeln. Wir arbeiten eng mit Spezialisten zusammen, um sicherzustellen, dass unsere Implantat-Systeme für Traumatologie/Orthopädie der oberen und unteren Extremitäten die Lebensqualität der Patienten verbessern. Die Sicherheit der Patienten steht bei uns immer an erster Stelle.

ES Trabajar con AGOMED significa trabajar con una empresa comprometida con la excelencia. Nuestros productos son diseñados en Alemania, biocompatibles y de vanguardia. Los miembros de nuestro equipo aportan, de manera colectiva, décadas de experiencia en tecnología médica.

Somos fiables. Somos flexibles. Somos innovadores.

AGOMED colabora con distribuidores, hospitales y doctores a nivel mundial con el objetivo de desarrollar soluciones rigurosas para problemas quirúrgicos complejos con tecnología punta. Cooperamos estrechamente con especialistas para garantizar que nuestros sistemas de implantes y instrumentos para Traumatología/Ortopedia de las extremidades superiores e inferiores mejoren la calidad de vida de los pacientes. La seguridad de los pacientes siempre es nuestra prioridad.

QUALITY IS PART OF OUR DNA



DIN EN ISO 13485

AGOMED MEANS QUALITY. WE EXCEED DIN EN ISO 13485 STANDARDS. ALL OUR PRODUCTS BEAR THE CE MARK.

EN



CE CERTIFICATE

AGOMED BEDEUTET QUALITÄT. WIR SIND ZERTIFIZIERT NACH DIN EN ISO 13485. ALLE UNSERE PRODUKTE TRAGEN DAS CE ZEICHEN.

DE

AGOMED ES SINÓNIMO DE CALIDAD. SUPERAMOS LOS ESTÁNDARES DIN EN ISO 13485. TODOS LOS PRODUCTOS LLEVAN EL SIGNO CE.

ES

INDEX

FEATURES & BENEFITS 2

ORDERING INFORMATION

TITANIUM PLATES 12

TITANIUM SCREWS 21

OPERATION TECHNIQUE

INDICATIONS 23

OPTION WITH SMALL GUIDE SLEEVES 26

USING CANNULATED SCREW 2.5 MM (OPTIONALITY) 34

OPTION WITH AIMING DEVICE 37

INSTRUMENTS FOR MULTIDIRECTIONAL FIXATION 42

INSTRUMENTS 45

SETS OVERVIEW 53

KING KONG

PLATING SYSTEM 2.5



LARGE SET S1464



SMALL SET S1486




PLATE FEATURES

ROUNDED EDGES &
A SMOOTH SURFACE

FOR SOFT TISSUE PROTECTION

2 DISTAL
SCREW ROWS

PROVIDE OPTIMAL
SUPPORT OF THE
ARTICULAR SURFACE

TRIPLE
CHAMFERED
PLATE EDGE

- VERY DISTAL PLATE POSITIONING POSSIBLE
- IMPROVED ANATOMICAL FIT ADAPTED TO THE VOLAR ASPECT OF THE DISTAL RADIUS
- MULTIDIRECTIONAL ANGULAR STABILITY OF $\pm 20^\circ$ IN ALL DIRECTIONS IN EACH SCREW HOLE

ANATOMICALLY PRE-CONTOURED PLATE DESIGN

LOW PLATE PROFILE

1.5 MM

DRILLING SLEEVES

ARE ALREADY PROVIDED INTO THE SCREW
HOLES TO ENSURE SAFE AND EASY DRILLING

COMPUTATIONAL FLUID DYNAMICS (CFD)

A PERFECTLY STREAMLINED SHAPE OF THE PLATE WAS OBTAINED THROUGH COMPUTATIONAL FLUID DYNAMICS WHICH IS USED IN AEROSPACE AND AUTOMOTIVE TECHNOLOGIES. THE DESIGN OF PLATES ENSURES THE FLOW LINES FOLLOW A SMOOTH TRAJECTORY, WITHOUT ABRUPT CHANGES IN DIRECTION

A PERFECTLY
STREAMLINED
SHAPE

ALL OUR PLATES HAVE A
MINIMAL EFFECT ON THE MUSCLES
AND DO NOT CAUSE PERCEPTIBLE
DISCOMFORT DURING MOVEMENT

PLATE FEATURES

HOOK PLATE

SPECIALLY MADE TO PREVENT OR TREAT NON-UNIONS AT THE ULNAR STYLOID



POINTED HOOKS GRIP

THE STYLOID PROCESS AND ACT AS REFERENCE POINT FOR PLATE APPLICATION

ULNAR HOOK PLATE

FOR CAPITAL AND SUBCAPITAL FRACTURES AND NON-UNIONS OF THE DISTAL ULNA



FOR DOUBLE PLATING TECHNIQUE

IN DORSALLY DISPLACED FRACTURES




RADIAL PLATE

LOW PLATE PROFILE
1.5 MM

- STABILIZES RADIAL COLUMN AND MAY BE USED FOR OTHER SMALL BONES
- EASILY CONTOURABLE TO PROVIDE THE DESIRED FIT

DORSAL PLATE



RADIAL PLATE

ANATOMICALLY PRE-CONTOURED



ULNAR PLATE



DORSAL PLATE

OBLONG HOLE & K-WIRE HOLES

AID IN PLATE POSITIONING

CLINICAL BENEFITS



WATERSHED LINE

ENHANCES PLACEMENT OF THE PLATE AND REDUCES THE RISK OF TENDON IRRITATIONS

FIRST DISTAL SCREW ROW

FOR SUPPORT OF THE CENTRAL ASPECT OF THE RADIOCARPAL JOINT

SECOND DISTAL SCREW ROW

PROVIDES STABILIZATION OF THE DORSAL RIM

MULTI-DIRECTIONAL SCREWS

LOCKING WITH A RANGE OF 20° OFFERS FRACTURE SPECIFIC PLATE FIXATION

SURFACE

TYPE II ANODIZATION (DOTIZE®) PROVIDES SMOOTH SURFACE, HIGHER FATIGUE RESISTANCE AND MINIMIZES THE RISK OF COLD-WELDING


K-WIRE HOLES

TO ASSIST WITH TEMPORARY PLATE FIXATION

- Ø 2.5 MM SCREWS CONSISTENT SCREW DIAMETER FOR ALL PLATES SIMPLIFIES THE INSTRUMENTS FOR EASE-OF-USE
- SCREWDRIVER TX7 FACILITATE SCREW INSERTION WITH SELF-HOLDING PROPERTIES AND SCREW REMOVAL WITH OPTIMUM TORQUE TRANSMISSION

ANATOMICAL PLATE DESIGN

SUPPORTS REDUCTION OF COMPLEX FRACTURES AND MINIMIZES THE NEED OF CONTOURING



BEST POSSIBLE EMBEDDING

IN SOFT TISSUE


ALL PLATES ARE AVAILABLE WITH & WITHOUT DRILLING SLEEVES



XXXXXXXX*
NUMBERS FOR PLATES COME WITH DRILLING SLEEVES

XXXXXXXX
NUMBERS FOR PLATES WITHOUT DRILLING SLEEVES

ADVANTAGES OF THE SCREWS



TAPERED SCREW HEAD

WITH DOUBLE THREAD PRODUCES A GRADUAL LOCKING EFFECT THAT PROMOTES A MORE SECURE SCREW INSERTION AND IMPROVED CONTOURING OF PLATE TO THE BONE



2X FASTER INSERTION WITH DOUBLE LEAD THREADS



ATRAUMATIC BLUNT TIP DESIGN

THE SCREW TIP CREATED AS TO BE MINIMALLY DAMAGING UPON INSERTION AS IN BICORTICAL FIXATION PROCEDURES




IMPROVED 6 HOLES SHAPE

IMPROVES THREADING CAPABILITIES OF THE LOCKING SCREW HEAD TO THE PLATE, EVEN IN OFF-AXIS CONFIGURATIONS



VARIABLE ANGLE LOCKING $\pm 20^\circ$



HEXALOBE DRIVE

THE SHAPE ENHANCES THE TORQUE THAT CAN BE APPLIED DURING INSERTION WHILE AT THE SAME TIME DECREASING THE CHANCES OF HEAD STRIPPING



ORDERING INFORMATION

KING TITANIUM PLATES / KING TITANPLATTEN / PLACAS DE TITANIO KING



NARROW / SCHMAL / ESTRECHA 



ADDITIONAL AIMING DEVICE / ZIEL GERÄT / BLOQUES GUÍA (ÁNGULO FIJO)





KONG TITANIUM PLATES / KONG TITANPLATTEN / PLACAS DE TITANIO KONG



BROAD / BREIT / ANCHA $\frac{1.5}{1.5}$ mm



4002600

4-holes, small, palmar, left
4-loch, klein, palmar, links
4-agujeros, pequeña, palmar, izquierda



4002601

4-holes, small, palmar, right
4-loch, klein, palmar, rechts
4-agujeros, pequeña, palmar, derecha



4002600*

4-holes, small, palmar, left
4-loch, klein, palmar, links
4-agujeros, pequeña, palmar, izquierda



4002601*

4-holes, small, palmar, right
4-loch, klein, palmar, rechts
4-agujeros, pequeña, palmar, derecha



4002602

5-holes, medium, palmar, left
5-loch, mittel, palmar, links
5-agujeros, mediana, palmar, izquierda



4002603

5-holes, medium, palmar, right
5-loch, mittel, palmar, rechts
5-agujeros, mediana, palmar, derecha



4002602*

5-holes, medium, palmar, left
5-loch, mittel, palmar, links
5-agujeros, mediana, palmar, izquierda



4002603*

5-holes, medium, palmar, right
5-loch, mittel, palmar, rechts
5-agujeros, mediana, palmar, derecha



4002604

6-holes, large, palmar, left
6-loch, groß, palmar, links
6-agujeros, grande, palmar, izquierda



4002605

6-holes, large, palmar, right
6-loch, groß, palmar, rechts
6-agujeros, grande, palmar, derecha



4002604*

6-holes, large, palmar, left
6-loch, groß, palmar, links
6-agujeros, grande, palmar, izquierda



4002605*

6-holes, large, palmar, right
6-loch, groß, palmar, rechts
6-agujeros, grande, palmar, derecha



ADDITIONAL AIMING DEVICE / ZIEL GERÄT / BLOQUES GUÍA (ÁNGULO FIJO)



1008056

Broad, left
Breit, links
Ancha, izquierda



1008057

Broad, right
Breit, rechts
Ancha, derecha



KING KONG TITANIUM PLATES / KING KONG TITANPLATTEN / PLACAS DE TITANIO KING KONG



LONG / LANGE / LARGA 2.0 mm



4002606

11-holes, XL, narrow, Palmar left
11-loch, XL, schmal, palmar, links
11-agujeros, XL, estrecha, palmar, izquierda

4002607

11-holes, XL, narrow, Palmar right
11-loch, XL, schmal, palmar, rechts
11-agujeros, XL, estrecha, palmar, derecha

4002606*

11-holes, XL, narrow, Palmar left
11-loch, XL, schmal, palmar, links
11-agujeros, XL, estrecha, palmar, izquierda

4002607*

11-holes, XL, narrow, Palmar right
11-loch, XL, schmal, palmar, rechts
11-agujeros, XL, estrecha, palmar, derecha



4002614

7-holes, L, narrow, palmar left
7-loch, L, schmal, palmar, links
7-agujeros, L, estrecha, palmar, izquierda

4002613

7-holes, L, narrow, palmar right
7-loch, L, schmal, palmar, rechts
7-agujeros, L, estrecha, palmar, derecha

4002614*

7-holes, L, narrow, palmar left
7-loch, L, schmal, palmar, links
7-agujeros, L, estrecha, palmar, izquierda

4002613*

7-holes, L, narrow, palmar right
7-loch, L, schmal, palmar, rechts
7-agujeros, L, estrecha, palmar, derecha

ADDITIONAL AIMING DEVICE / ZIEL GERÄT / BLOQUES GUÍA (ÁNGULO FIJO)



1008058

L / XL, narrow, left
L / XL, schmal, links
Estrecha L / XL, izquierda



1008059

L / XL, narrow, right
L / XL, schmal, rechts
Estrecha L / XL, derecha



**KING KONG TITANIUM PLATES / KING KONG TITANPLATTEN /
PLACAS DE TITANIO KING KONG**



HOOK $\frac{1.5 \text{ mm}}$



4002612
6-holes, hook
6-loch, hacken
6-agujeros,
gancho

4002612*
6-holes, hook
6-loch, hacken
6-agujeros,
gancho

STRAIGHT /
GERADE / RECTA $\frac{1.5 \text{ mm}}$

L FORM $\frac{1.5 \text{ mm}}$

T FORM $\frac{1.5 \text{ mm}}$



4002610
6-holes
6-loch
6-agujeros



4002610*
6-holes
6-loch
6-agujeros



4002609
2/4-holes, right
2/4-loch, rechts
2/4-agujeros,
derecha



4002608
2/4-holes, left
2/4-loch, links
2/4-agujeros,
izquierda



4002609*
2/4-holes, right
2/4-loch, rechts
2/4-agujeros,
derecha



4002608*
2/4-holes, left
2/4-loch, links
2/4-agujeros,
izquierda



4002632
7-holes
7-loch
7-agujeros



4002611
8-holes
8-loch
8-agujeros



4002633
10-holes
10-loch
10-agujeros



4002632*
7-holes
7-loch
7-agujeros



4002611*
8-holes
8-loch
8-agujeros



4002633*
10-holes
10-loch
10-agujeros



SAMPLE INDICATION:
• Fractures and Osteotomies of the distal ulna



ATRAUMATIC FRAME TITANIUM PLATES / ATRAUMATISCHE RAHMENPLATTE / MARCO ATRAUMÁTICO DORSAL DE TITANIO



SAMPLE INDICATION:
• Multiple fractures on the distal radius with dorsal incision

1.5 mm



4002615
14-holes, size 1, right
14-loch, groesse 1, rechts
14-agujeros, tamaño 1, derecha



4002616
14-holes, size 1, left
14-loch, groesse 1, links
14-agujeros, tamaño 1, izquierda



4002617
14-holes, size 2, right
14-loch, groesse 2, rechts
14-agujeros, tamaño 2, derecha



4002618
14-holes, size 2, left
14-loch, groesse 2, links
14-agujeros, tamaño 2, izquierda


TITANIUM SCREWS 2.5 / TITANSCHRAUBEN 2.5 / TORNILLOS DE TITANIO 2.5



- Double Thread
- 200 % Faster implantation

Length / Länge / Longitud
Multidirectional, Locking, TX 7 / Multidirectional, Winkelstabil, TX 7 / Multidirectional, Cabeza roscada, TX 7
Self-tapping, Standard, TX 7 / Selbstschneidend, Standard, TX 7 / Autoroscantes, Estandar TX 7

L	Ø 2.5 mm	Ø 2.5 mm	Pack/Pack/Paq.
10 mm	5028010	5028110	5
11 mm	5028011	5028111	5
12 mm	5028012	5028112	5
13 mm	5028013	5028113	5
14 mm	5028014	5028114	5
15 mm	5028015	5028115	5
16 mm	5028016	5028116	5
17 mm	5028017	5028117	5
18 mm	5028018	5028118	5
19 mm	5028019	5028119	5
20 mm	5028020	5028120	5
21 mm	5028021	5028121	5
22 mm	5028022	5028122	5
23 mm	5028023	5028123	5
24 mm	5028024	5028124	5
25 mm	5028025	5028125	5
26 mm	5028026	5028126	5
27 mm	5028027	5028127	5
28 mm	5028028	5028128	5
29 mm	5028029	5028129	5
30 mm	5028030	5028130	5



Atraumatic frame titanium plates are available with drilling sleeves (4002615*, 4002616*, 4002617*, 4002618*)




**CANNULATED TITANIUM SCREWS 2.5 / KANÜLIERTE TITANSCHRAUBEN 2.5 /
TORNILLOS DE TITANIO CANULADO 2.5**



L	Ø 2.5 mm	Pack/Pack/Paq.
10 mm	5028210	5
11 mm	5028211	5
12 mm	5028212	5
13 mm	5028213	5
14 mm	5028214	5
15 mm	5028215	5
16 mm	5028216	5
17 mm	5028217	5
18 mm	5028218	5
19 mm	5028219	5
20 mm	5028220	5
21 mm	5028221	5
22 mm	5028222	5
23 mm	5028223	5
24 mm	5028224	5
25 mm	5028225	5
26 mm	5028226	5
27 mm	5028227	5
28 mm	5028228	5
29 mm	5028229	5
30 mm	5028230	5

-  Drill-ø 2.0 mm
-  K-wire-ø 1.1 mm





OPERATION TECHNIQUE

INDICATIONS

- Fixation of complex intra- and extra-articular fractures and osteotomies of the distal radius
- Fractures and osteotomies of the distal ulna
- Correction osteotomies



CONTRAINDICATIONS

- Infection or inflammation (local or systemic)
- Allergies to the implant material
- Acute or chronic osteomyelitis at or close to the surgical field
- Insufficient soft tissue coverage
- Severe soft tissue swelling compromising normal wound healing
- Unacceptably high anesthesia risk
- Fractures in children and adolescents with epiphyseal plates that are not yet ossified



SAMPLE INDICATION:

- Fixation of complex intra- and extra-articular fractures and osteotomies of the distal radius



OPERATION TECHNIQUE OPTION WITH SMALL GUIDE SLEEVES

Preoperative planning

- In addition to taking standard X-rays in A/P and sagittal planes in neutral position of the wrist, a high-resolution computer tomography should be conducted for the further clarification of intra-articular fractures.

Patient Positioning


- Place the patient in the supine position with the hand and arm on a hand table, preferably radiolucent for fluoroscopic imaging. The elbow should be fully extended and in full supination and under tourniquet control.
- A cloth roll can be placed underneath the wrist as a reduction aid.

Incision

- Make a longitudinal incision slightly radial to the flexor carpi radialis tendon (FCR).

Soft tissue preparation / Exposure of the radius

- Retract the FCR tendon toward the ulna while protecting the median nerve. Dissect through the floor of the FCR sheath, protecting the radial artery radially, exposing the pronator quadratus. Release the pronator quadratus radially to ulnarly by sharply incising over the watershed line and proximally on the lateral edge of the radius.
- Elevate the PQ to expose the volar surface of the radius. The origin of the flexor pollicis longus (FPL) muscle can be partially released if needed for added exposure.



Precaution:

Leave the volar wrist capsule intact to avoid devascularization of the fracture fragments and destabilization of the volar wrist ligaments.

Fracture Reduction

Before Plating:

- Reduce the fracture using the preferred reduction technique under fluoroscopic imaging. The reduction method will be fracture specific.
- Insert K-wires to hold the reduction.
- If needed, K-wires can be inserted parallel to the subchondral bone to reconstruct the articular surface. Ensure that K-wires are placed in locations which will not impede proper plate placement and fixation.
- Confirmed under fluoroscopy in AP and lateral views.

STEP 1



Insertion of the plate

- Select the plate according to fracture type and radial anatomy. To conform different fracture types, bone dimensions and bone qualities, a choice of narrow and broad volar plates is offered to the surgeon.

- The plate is designed to sit along the distal aspect of the radius to support the volar articular fracture fragments and just proximal to the watershed line to minimize tendon irritation.

The plate can be adjusted proximal or distal through the use of the oblong hole.

- Particular emphasis must be placed on restoration of the correct length correlation between radius and ulna, elimination of any steps in the radial joint surface, and compensation of any radial offset or malrotation of the distal fragment since otherwise this might restrict the postoperative range of motion.



Note:

If necessary, insert K-wires $\phi 1.4$ through the plate to preliminarily fixate it. K-wire holes are in the distal and the shaft part of the plate. They may help to preliminarily fixate the result of reduction, depending on the fracture pattern.



Plate Bending

- If required, the frame plates, the small fragment plates, the hook plates and the distal ulna plates can be bent with the plate bending pliers.
- The labeled side of the plate must always face upwards when inserting the plate into the bending pliers.
- While bending, the plate must always be held at two adjacent holes to prevent contour deformation of the intermediate plate hole.



Note:

Do not bend the plate by more than 30°. Bending the plate further may deform the plate holes and may cause the plate to break postoperatively.

STEP 2



• Drilling the core hole

After fracture reduction, the first borehole is made into the slotted hole of the shaft using the drill guide and the twist drill $\phi 2.0$.



STEP 3



• Determination of the screw length

Measure screw length using the depth gauge.

• Insertion of the first shaft screw

Insert a non-locking screw $\phi 2.5$ into the oblong shaft hole, using the self-retaining TX7 screwdriver.

Note:

Positioning slot only accepts non-locking screw $\phi 2.5$.

- Confirm plate position by fluoroscopy in A/P and lateral views. Particularly consider correct placement of the distal part of the plate and the alignment of the plate shaft. Adjust improper plate position if required.
- If necessary, the result has to be corrected and the plate displaced in longitudinal and/or lateral direction. The screw has to be loosened for this purpose.



Note:

To avoid contact with flexor tendons, the plate must be applied just proximal to and below the watershed line.



STEP 4



- **Insertion of another shaft screw**
In order to be able to absorb optimally the forces in the shaft region during reduction, it is advisable to insert another shaft screw, a locking screw, prior to the reduction, ensuring that the plate is positioned correctly.
- Place a locking screw $\phi 2.5$ into the hole using the drill guide and a drill $\phi 2.0$.




Note:

The screwable drill-guide (1008065) can be locked with a clockwise turn in the holes of the plate. It thus performs all of the functions of a drill guide without the need to be held. (see step 7)

- Locking screws may be inserted fixed-angled (0°) or with variable angle. An off-center angle exceeding 20° is not recommended as it may prevent the screw from locking correctly into the plate and entail screw loosening.
- Except for the oblong hole, all other plate holes can accommodate locking as well as non-locking screws with an angulation of up to 20° off-axis.
- The tips of the screws should not protrude beyond the dorsal cortex since otherwise this might injure the extensor tendons. As the screws are supplied in 2 mm increments, the ideal screw length corresponding to the measurement result may not be available. Therefore, select the screw with the next shorter length.

Note:

When combining locking and non-locking screws always insert the non-locking screw first, in order to pull the plate to the bone.



- To remove the screws from the implant container, insert the appropriate screwdriver perpendicularly into the screw head of the desired screw.
- Vertically extract the screw from the compartment.
- Check the screw length at the scale of the measuring module.

STEP 5



- The fracture reduction can be fixed with temporary K-wires through the distal holes.
- The K-wires can be positioned in such a way that the position of the plate to the distal radioulnar joint as well as to the radiocarpal joint can be checked simultaneously.

Note:

The K-wires can be positioned through the small drill sleeves, by using the sleeve reduction guide (1008069)





STEP 6



- **Distal plate fixation**
Drill the first core holes (Ø 2.0 mm) for the locking screws (Ø 2.5 mm) and measure the screw length.
- Facilitate fixed angle drilling by using blue small drill sleeves, mounted onto the plate.
- Drill to the desired depth with a drill ø2.0.
- Check with fluoroscopy.

Note:

The instruments are offered for **multidirection screw fixation**: a threaded funnel-shaped version (1008041) and drill-guides (1008040 & 1008066). They must be used without blue small drill sleeves (see page 42).



- Remove the small drill sleeve, determine screw length with the depth gauge.
- The ideal screw length for the assigned drill hole can be read on the scale of the depth gauge.

Note:


The depth gauge caliper has a hooked tip that is either inserted to the bottom of the hole or is used to catch the far cortex of the bone.



- Insert an appropriate length locking screw using the TX7 screwdriver stem adjusting.

Note:

Do not overtighten the screw.



- Using the 2.0 mm bit, drilling by using blue small drill sleeves in series all remaining available holes.
- Remove the blue small drill sleeves, determine screw length with the depth gauge.
- Insert an appropriate length locking screw, using the TX7 screwdriver stem adjusting until final tightening.

Note:

Each hole should be prepared sequentially.



- Insert further screws according to fracture pattern. It is recommended to fully accommodate the distal screw row and place screws in the second row.
- Once all distal screws have been inserted, remove the K-wires.




STEP 7



- The most proximal plate hole should be preferred when securing the plate shaft with not less than 2 locking screws.
- Insert the screwable drill-guide (1008065) and drill using 2.0mm drill.

Note:

Can be locked with a clockwise turn in the holes of the plate.



- Finally, confirm all screw heads have locked flush with the plate surface.
- Drill all remaining available holes and insert locking screws.

USING CANNULATED SCREW 2.5 MM (OPTIONALITY)



- Insert the sleeve reduction guide (1008069) into the desired hole with a small drill sleeve. The sleeve reduction guide will center the 1.1 mm K-wire during insertion.
- Confirm K-wire placement under fluoroscopy.

Note:

Use 1008069 together with the tool 1008009



- Drill using the Cannulated Drill (1019106).
- Remove the small drill sleeve, determine screw length with the depth gauge.
- Insert the cannulated screw 2.5 mm over K-wire using the TX7 cannulated screwdriver (1001911).

Note:

Cannulated screwdriver is not for final tightening of the locking screw



- Once screw engages the plate, remove the K-wire and finish final tightening of the cannulated screw using the solid core TX7 screwdriver stem (1001910).



- **Cannulated Variable Angle Screw Placement**
Cannulated Locking Screws 2.5 mm can be inserted at $\pm 20^\circ$.





STEP 8



- Check the final result under fluoroscopy in A/P and lateral views.
- Exchange dorsally protruding screws for shorter ones and correct mal-positioned screws, if necessary.

Note:


The screws in the first row should slightly be tilted proximally; by contrast, the screws in the second row should be tilted distally. The subchondral screw arrangement according to the array principle provides optimum support for both the central region and the dorsal edge of the radius.

Wound closure

- The wound is closed in layers. Following the skin closure, a final x-ray image is taken.

Postoperative treatment

- After completion of the surgery, a forearm splint is applied, which allows active finger movement.
- In general, the follow-up treatment is individually defined by the surgeon, and depends on fracture type and bone quality.



OPERATION TECHNIQUE OPTION WITH AIMING DEVICE

STEP 1



- The volar radial plate can be applied as a buttress plate in flexion fractures, but may also address extension fractures.
- Fracture reduction is gradually performed with subsequent retention by locking plate fixation.
- Particular emphasis must be placed on restoration of the correct length correlation between radius and ulna, elimination of any steps in the radial joint surface, and compensation of any radial offset or malrotation of the distal fragment since otherwise this might restrict the postoperative range of motion. Hence, careful fracture reduction with X-ray control is crucial. K-wires or reduction forceps may help to preliminarily fixate the result of reduction, depending on the fracture pattern.



Note:

With the wrist in hyper-flexion lateral fluoroscopy will aid in obtaining a precise sub-chondral position of screws in the distal fragment to the extent possible. This type of fixation ensures the best placement for the screws and helps to avoid secondary dislocation.

- To conform different fracture types, bone dimensions and bone qualities, a choice of narrow and broad volar plates is offered to the surgeon.





STEP 2



- Place the plate proximal to the distal edge of the radius. Following the watershed line, the distal part of the plate should be lower on the radial side.
- If necessary, insert K-wires $\phi 1.4$ through the plate to preliminarily fixate it. K-wire holes are in the distal and the shaft part of the plate.

STEP 3



- **Drilling the core hole**
After fracture reduction, the first borehole is made into the slotted hole of the shaft using the drill guide (1008066) and the twist drill $\phi 2.0$.

STEP 4



- **Determination of the screw length**
Measure screw length using the depth gauge.
- **Insertion of the first shaft screw**
Insert a non-locking screw $\phi 2.5$ into the oblong shaft hole, using the self-retaining TX7 screwdriver.

Note:

Positioning slot only accepts non-locking screw $\phi 2.5$.



- Confirm plate position by fluoroscopy in A/P and lateral views. Particularly consider correct placement of the distal part of the plate and the alignment of the plate shaft. Adjust improper plate position if required.
- If necessary, the result has to be corrected and the plate displaced in longitudinal and/or lateral direction. The screw has to be loosened for this purpose.

Note:


To avoid contact with flexor tendons, the plate must be applied just proximal to and below the watershed line.



STEP 5



- **Insertion of screws**
Locking screws may be inserted fixed-angled (0°) or with variable angle. An off-center angle exceeding 20° is not recommended as it may prevent the screw from locking correctly into the plate and entail screw loosening.



- Facilitate fixed angle drilling by using an appropriate targeting device mounted onto the plate with the fixing screw.
- Engage a fixed angle drill guide (1008065) in the selected distal plate hole and drill to the desired depth with a drill $\phi 2.0$.



Note:
The aiming devices are designed for seating fixed-angled locking screws only.

STEP 6



- Remove the drill guide, determine screw length with the depth gauge (1006002).
- Repeat these steps for all other distal holes.

IMPORTANT:

While measuring the lengths with a depth guide, the height of the aiming device must be considered, for that needed from obtained height of the depth gauge subtract **6 mm**.



- Insert an appropriate length locking screw using the TX7 screwdriver stem adjusting.
- Once all distal screws have been inserted, remove the aiming device.
- The most proximal plate hole should be preferred when securing the plate shaft with not less than 2 locking screws.
- Finally, confirm all screw heads have locked flush with the plate surface.
- Drill all remaining available holes and insert locking screws.







INSTRUMENTS FOR MULTIDIRECTIONAL SCREW FIXATION



- The instruments are offered for variable angled drilling: a threaded funnel-shaped version (1008041) and drill-guides (1008066 & 1008040).



- The drill-guides 1008066 & 1008040 are intended for stabilization and keeping the drill at a certain angle while drilling holes.
- Holding the drill-guide 1008009 direct the drill to the required hole at the appropriate angle.



Note:

The drill-guides 1008066 & 1008040 allow for a multidirectional angulation of $\pm 20^\circ$, so that fixed-angle locking is always ensured.



- Use one of the variable angle drill-guides in a locking plate hole and drill with a drill $\phi 2.0$ in a chosen angle.
- Use the depth gauge to confirm screw length.
- Insert selected screw using the self-retaining TX7 screwdriver stem adjusting until final tightening.



INSTRUMENTS FOR MULTIDIRECTIONAL K-WIRE FIXATION



- The K-wire guide tools 1008068 & 1008050 are intended for stabilization and keeping k-wire at a certain angle during its insertion. ($\pm 20^\circ$)
- Holding the tool 1008009, direct the k-wire at the appropriate angle to the required hole.
- Use together with the tool 1008009



for K-wire $\phi 1.1\text{mm}$



STEP 7



- Check the final result under fluoroscopy in A/P and lateral views.
- Exchange dorsally protruding screws for shorter ones and correct mal-positioned screws, if necessary.

Note:

The screws in the first row should slightly be tilted proximally; by contrast, the screws in the second row should be tilted distally. The subchondral screw arrangement according to the array principle provides optimum support for both the central region and the dorsal edge of the radius.

Wound closure

- The wound is closed in layers. Following the skin closure, a final x-ray image is taken.

Postoperative treatment

- After completion of the surgery, a forearm splint is applied, which allows active finger movement.
- In general, the follow-up treatment is individually defined by the surgeon, and depends on fracture type and bone quality.



INSTRUMENTS

1000321

Agopaqx4 tray for implant modules w/4 lids, King Kong system, w/o implant modules

Agopaq halte-tray f. implantate, mit 4 deckeln, King Kong system, ohne implantatmodule

Agopaqx4 soporte p/modules con implantes c/4 tapas, King Kong sistema, sin modulos p/implantes



1301037

Agopaqx2 tray for implant modules w/ 2 lids, w/o implant modules

Agopaqx2 halte-tray f. implantate, mit 2 deckeln, ohne implantatmodule

Agopaqx2 soporte p/modules con implantes c/2 tapas, sin modulos p/implantes




1000376

King Kong sterilisation tray for instruments, w/o instruments

King Kong sterilisiersiebschale f. instrumente, ohne instrumente

King Kong bandeja para instrumental, sin instrumental





1000438

Agopaqx2 instrument tray f.
small system 2.5, w/o instruments
*Agopaqx2-instrumenten-siebschale
f. klein system 2.5, ohne instrumente
Agopaqx2 caja para instrumental
p/sistema pequeña 2.5, sin instrumental*



1000380

Agopaqx4 plate-tray King Kong
system 2.5, Broad, w/o implants
*Agopaqx4 platten-modul King Kong
system 2.5, breit, ohne implantate
Agopaqx4 soporte p/placas King Kong
sistema 2.5, Ancha, sin implantes*



1000381

Agopaqx4 plate-tray King Kong
system 2.5, w/o implants
*Agopaqx4 platten-modul King Kong
system 2.5, ohne implantate
Agopaqx4 soporte p/placas King Kong
sistema 2.5, sin implantes*



1000384

Agopaqx4 screw tray King Kong
system 2.5, w/o implants
*Agopaqx4 schrauben modul
King Kong system 2.5, ohne implantate
Agopaqx4 soporte p/tornillos King Kong
sistema 2.5, sin implantes*



1000379

Agopaqx4 plate-tray King Kong
system 2.5, Narrow, w/o implants
*Agopaqx4 platten-modul King Kong
system 2.5, schmal, ohne implantate
Agopaqx4 soporte p/placas King Kong
sistema 2.5, Estrecha, sin implantes*



1000382

Agopaqx4 plate-tray King Kong
system 2.5, Long XL, w/o implants
*Agopaqx4 platten-modul King Kong
system 2.5, Lange XL, ohne implantate
Agopaqx4 soporte p/placas King Kong
sistema 2.5, Larga XL, sin implantes*



1000322

Agopaqx4 tray King Kong system 2.5,
Atraumatic frame plate, w/o implants
*Agopaqx4 platten-modul King Kong system 2.5,
atraumatische rahmenplatte, ohne implantate
Agopaqx4 soporte p/placas King Kong sistema 2.5,
placa marco atraumática, sin implantes*



1000383

Agopaqx4 tray King Kong system 2.5,
Aiming Device, w/o implants
*Agopaqx4 modul King Kong system 2.5,
zielgerät, ohne implantate
Agopaqx4 soporte p/dispositivo de puntería,
King Kong sistema 2.5, sin implantes*





1019106

Twist drill 2,0 x 95 mm, cannulated 1,2 mm,
AO-shaft, stop 26 mm

*Spiralbohrer 2,0 x 95 mm,
kanueltiert 1,2 mm AO, nl: 26 mm
Broca 2,0 x 95 mm, canulada 1,2 mm
vastago AO, stop 26 mm*



1201008

Twist drill 2,0 x 103 mm,
AO-shaft, grey code

*Spiralbohrer 2,0 x 103 mm,
59 mm NL, AO-ansatz, farbcode grau
Broca 2,0 x 103 mm,
vastago AO, codido gris*



1019103

Twist drill 2,0 x 103 mm,
AO-shaft, grey code

*Spiralbohrer 2,0 x 95 mm,
kanueltiert 1,2 mm rundschaft, nl: 26 mm
Broca 2,0 x 95 mm, canulada 1,2 mm
redonda, stop 26 mm*



i OPTIONAL

1201029

Twist drill 1,8 x 95 mm,
stop 29 mm, AO, pink code

*Spiralbohrer 1,8 x 95 mm,
stop 29 mm, AO, farbcode pink
Broca 1,8 x 95 mm,
stop 29 mm, AO, codigo rosa*



1001910

Screwdriver blade, double TX 7
non-cannulated, 90 mm, AO-shaft

*Schraubendreherklinge, doppelt TX 7
nicht kanueltiert, 90 mm, AO-ansatz
Hoja intercambiable, doble TX 7
no canulada, 90 mm, vastago AO*



1001911

Screwdriver blade, double TX 7
cannulated, 90 mm, AO-shaft

*Schraubendreherklinge, doppelt TX 7
kanueltiert, 90 mm, AO-ansatz
Hoja intercambiable, doble TX 7
canulada, 90 mm, vastago AO*



1004036

Silicone screwdriverhandle,
black cannulated, 12 cm, AO-shaft

*Silikon-schraubendrehergriff schwarz
kanueltiert gross, 12 cm, AO-ansatz
Mango destornillador silicona canulado,
grande, 12 cm, negro, vastago-AO*



1004049

Silicone screwdriverhandle,
black cannulated, 15 cm, AO-shaft

*Silikon-schraubendrehergriff schwarz
kanueltiert gross, 15 cm, AO-ansatz
Mango destornillador silicona canulado,
grande, 15 cm, negro, vastago-AO*



1002527

Plate- /screwholding forceps
angled, 15cm/6"

*Titan platten-/schraubenhaltepinzette
15cm/6" gewinkelt
Pinza de titanio p/placas y tornillos
15cm/6" angulada*



1006024

Agomed depth gauge 40 mm,
round f. 2,0 - 4,5 mm screws

*Agomed tiefenlehre 40 mm, rund
f. 2,0 - 4,5 mm schrauben, taster 1,6 mm
Agomed medidor profundidad 40 mm
p/tornillos 2,0 - 4,5 mm*



1008009

Handle for drill-guide,
only hand-systems

*Bohrbuechsengriff allein
f. hand-systeme
Mango solo p/guia de
broca sistemas mano*



1005620

Flat bending plier, 13,5cm/5 1/2"

*Flachbiegezange 13,5cm/5 1/2"
Pinza doblar placas, bocas planas,
13,5cm/5 1/2"*





1101340

Olive kirschner-wires, thread 10 mm,
trocar/round, 60 x 1.1 mm, pack/1
*Olive kirschner draht, gewinde 10 mm,
trokar/rund, 60 x 1,1 mm, pack/1*
*Olive broca filiforme kirschner, rosca 10 mm,
trocar/redondo, 60 x 1,1 mm, pack/1*



1101023

Kirschner wire 100 x 1,1 mm
trocar/round, pack/6
*Kirschner draht 100 x 1,1 mm
trokar/rund, pack/6*
*Kirschner agujas 100 x 1,1 mm
trocar/redondo, paq./6*



1101063

Werber CBS-K-wire 150 x 1,4 mm
trocar / round, pack/6
*Werber CBS-Kirschner draht
150 x 1,4 mm, trokar/rund, pack/6*
*Werber CBS-Kirschner agujas
150 x 1,4 mm, trocar/redondo, paq./6*



1104009

Agomed k-wire dispenser 15 cm/6"
for k-wire 0,8 - 1,6 mm diameter
*Agomed kirschnerdrahtspender 15 cm/6"
f. drahte 0,8 - 1,6 mm durchmesser*
*Agomed dispensador agujas kirschner
15 cm/6", p/agujas de 0,8 - 1,6 mm*



1008069

Sleeve reduction guide,
K-wire $\varnothing$ 1.1, length 35 mm
*Sleeve reduction guide
K-draht $\varnothing$ 1.1, lnge 35 mm*
*Gu de reducci3n de mangas,
kirschner agujas $\varnothing$ 1.1, longitud 35 mm*



1008041

Drill guide with thread,
drill $\varnothing$ 2.0, 0° to 10°
*Bohrhlse mit gewinde,
bohren $\varnothing$ 2.0, 0°-10°*
*Gu con rosca,
broca $\varnothing$ 2.0, 0°-10°*



1008065

Drill-guide only, screwable,
length 45 mm, drill $\varnothing$ 2.0
*Bohrlehren-aufsatz allein,
schraubbar 45 mm, bohren $\varnothing$ 2.0*
*Gu con rosca, p/broca
longitud 45 mm, broca $\varnothing$ 2.0*



1008051

K-wire guide only, screwable,
K-wire $\varnothing$ 1.1, length 45 mm
*Bohrlehren-aufsatz allein,
schraubbar K-draht $\varnothing$ 1.1, lnge 45 mm*
*Gu con rosca, p/kirschner
agujas $\varnothing$ 1.1, longitud 45 mm*



1008066

Drill-guide only, variable,
length 35 mm, drill $\varnothing$ 2.0
*Bohrlehren-aufsatz allein,
variabel lnge 35 mm, bohren $\varnothing$ 2.0*
*Gu broca p/tornillos,
variable, longitud 35 mm, broca $\varnothing$ 2.0*



1008068

K-wire guide only, variable,
K-wire $\varnothing$ 1.1, length 35 mm
*Bohrlehren-aufsatz allein,
variabel K-draht $\varnothing$ 1.1, lnge 35 mm*
*Gu variable, p/kirschner agujas $\varnothing$ 1.1,
longitud 35 mm*



1008040

Drill-guide only, variable,
length 35 mm, drill $\varnothing$ 2.0
*Bohrlehren-aufsatz allein,
variabel lnge 35 mm, bohren $\varnothing$ 2.0*
*Gu broca p/tornillos,
variable, longitud 35 mm, broca $\varnothing$ 2.0*



1008050

K-wire guide only, variable,
K-wire $\varnothing$ 1.1, length 35 mm
*Bohrlehren-aufsatz allein,
variabel K-draht $\varnothing$ 1.1, lnge 35 mm*
*Gu variable, p/kirschner agujas $\varnothing$ 1.1,
longitud 35 mm*





1008054

Aiming device Radius 2.5,
narrow, left
*Ziel gerät Radius 2.5,
schmal, links*
*Guia para Radio 2.5,
estrecha, izquierda*



1008055

Aiming device Radius 2.5,
narrow, right
*Ziel gerät Radius 2.5,
schmal, rechts*
*Guia para Radio 2.5,
estrecha, derecha*



1008056

Aiming device Radius 2.5,
broad, left
*Ziel gerät Radius 2.5,
breit, links*
*Guia para Radio 2.5,
ancha, izquierda*



1008057

Aiming device Radius 2.5,
broad, right
*Ziel gerät Radius 2.5,
breit, rechts*
*Guia para Radio 2.5,
ancha, derecha*



1008058

Aiming device Radius 2.5,
Long L / XL, left
*Ziel gerät Radius 2.5,
lange L / XL schmal, links*
*Guia para Radio 2.5, larga
estrecha L / XL, izquierda*



1008059

Aiming device Radius 2.5,
Long L / XL, right
*Ziel gerät Radius 2.5,
lange L / XL schmal, rechts*
*Guia para Radio 2.5, larga
estrecha L / XL, derecha*



S1464 - KING KONG SYSTEM 2.5 SET



Pos.	Cat.No.	Description	Qty.	Picture page
Containers				
1	1000321	AGOPAQX4 TRAY FOR IMPLANT MODULES W/4 LIDS, KING KONG SYSTEM, W/O IMPLANT MODULES	1	45
2	1000322	AGOPAQX4 PLATE-TRAY KING KONG SYSTEM 2.5, ATTRAUMATIC FRAME PLATE, W/O IMPLANTS	1	47
3	1000376	KING KONG STERILISATION TRAY FOR INSTRUMENT, W/O INSTRUMENTS	1	45
4	1000379	AGOPAQX4 PLATE-TRAY SYSTEM 2.5, NARROW, W/O IMPLANTS	1	47
5	1000380	AGOPAQX4 PLATE-TRAY KING KONG SYSTEM 2.5, BROAD, W/O IMPLANTS	1	46
6	1000381	AGOPAQX4 PLATE-TRAY SYSTEM 2.5, W/O IMPLANTS	1	46
7	1000382	AGOPAQX4 PLATE-TRAY KING KONG SYSTEM 2.5, LONG XL, W/O IMPLANTS	1	47
8	1000383	AGOPAQX4 TRAY KING KONG SYSTEM 2.5, AIMING DEVICE, W/O IMPLANTS	1	47
9	1000384	AGOPAQX4 SCREW TRAY SYSTEM 2.5, W/O IMPLANTS	1	46
Instruments				
10	1001910	SCREWDRIVER BLADE, DOUBLE TX 7 NON-CANNULATED, 90 MM, AO-SHAFT	1	48
11	1001911	SCREWDRIVER BLADE, DOUBLE TX 7 CANNULATED, 90 MM, AO-SHAFT	1	48
12	1004049	SILICONE SCREWDRIVERHANDLE BLACK, CANNULATED, 15 CM, AO-SHAFT	1	49
13	1004036	SILICONE SCREWDRIVERHANDLE BLACK CANNULATED, 12 CM, AO-SHAFT	1	49
14	1005620	FLAT BENDING PLIER, 13,5CM/5 1/2"	2	49
15	1006024	AGOMED DEPTH GAUGE 40 MM, ROUND F. 2,0 - 4,5 MM SCREWS	1	49
16	1008009	HANDLE FOR DRILL-GUIDE, ONLY HAND-SYSTEMS	1	49

Pos.	Cat.No.	Description	Qty.	Picture page
17	1201008	TWIST DRILL 2,0 X 103 MM, AO-SHAFT, GREY CODE	1	48
18	1019103	TWIST DRILL 2,0 X 95 MM, CANNULATED 1,2 MM ROUND SHAFT, STOP 26 MM	1	48
19	1008041	DRILL GUIDE ONLY, SCREWABLE, DRILL Ø2.0, 0° TO 10°, SYSTEM 2.5	1	50
20	1002527	PLATE- /SCREWHOLDING FORCEPS ANGLED, 15CM/6"	1	49
21	1008040	DRILL-GUIDE ONLY, VARIABLE, LENGTH 35 MM, DRILL Ø 2.0, ROUND TIP	1	51
22	1008050	K-WIRE GUIDE ONLY, VARIABLE, K-WIRE Ø 1.1, LENGTH 35 MM, ROUND TIP	1	51
23	1008051	K-WIRE GUIDE ONLY, SCREWABLE, K-WIRE Ø 1.1, LENGTH 45 MM	1	51
24	1008065	DRILL-GUIDE ONLY, SCREWABLE, LENGTH 45 MM, DRILL Ø 2.0	1	51
25	1008066	DRILL-GUIDE ONLY, VARIABLE, LENGTH 35 MM, DRILL Ø 2.0	1	51
26	1008068	K-WIRE GUIDE ONLY, VARIABLE, K-WIRE Ø 1.1, LENGTH 35 MM	1	51
27	1008069	SLEEVE REDUCTION GUIDE, 35 MM LENGTH, K-WIRE Ø1,1, SYSTEM 2.5	1	50
28	1008054	AIMING DEVICE RADIUS 2.5, NARROW, LEFT	1	52
29	1008055	AIMING DEVICE RADIUS 2.5, NARROW, RIGHT	1	52
30	1008056	AIMING DEVICE RADIUS 2.5, BROAD, LEFT	1	52
31	1008057	AIMING DEVICE RADIUS 2.5, BROAD, RIGHT	1	52
32	1008058	AIMING DEVICE RADIUS 2.5, LONG L / XL NARROW, LEFT	1	52
33	1008059	AIMING DEVICE RADIUS 2.5, LONG L / XL NARROW, RIGHT	1	52
34	1019106	TWIST DRILL 2,0 X 95 MM, CANNULATED 1,2 MM, AO-SHAFT, STOP 26 MM	1	48



Pos.	Cat.No.	Description	Qty.	Picture page
35	1101340	OLIVE KIRSCHNER-WIRES, THREAD 10 MM, TROCAR/ROUND, 60 X 1.1 MM, PACK/1	3	50
36	1101023	KIRSCHNER WIRE 100 X 1,1 MM TROCAR/ROUND, PACK/6	1	50
37	1101063	WERBER CBS-K-WIRE 150 X 1,4 MM TROCAR / ROUND, PACK/6	1	50
38	1104009	AGOMED K-WIRE DISPENSER 15 CM/6" FOR K-WIRE 0,8 - 1,6 MM DIAMETER	2	50
Plates				
39	4002594	TITANIUM PLATE, 4-HOLES, NARROW SMALL, PALMAR, LEFT, SYSTEM 2.5	1	12
40	4002595	TITANIUM PLATE, 4-HOLES, NARROW SMALL, PALMAR, RIGHT, SYSTEM 2.5	1	12
41	4002596	TITANIUM PLATE, 5-HOLES, NARROW MEDIUM, PALMAR, LEFT, SYSTEM 2.5	1	13
42	4002597	TITANIUM PLATE, 5-HOLES, NARROW MEDIUM, PALMAR, RIGHT, SYSTEM 2.5	1	13
43	4002598	TITANIUM PLATE, 6-HOLES, NARROW LARGE, PALMAR, LEFT, SYSTEM 2.5	1	13
44	4002599	TITANIUM PLATE, 6-HOLES, NARROW LARGE, PALMAR, RIGHT, SYSTEM 2.5	1	13
45	4002600	TITANIUM PLATE, 4 - HOLES, BROAD SMALL, PALMAR, LEFT, SYSTEM 2.5	1	14
46	4002601	TITANIUM PLATE, 4 - HOLES, BROAD SMALL, PALMAR, RIGHT, SYSTEM 2.5	1	14
47	4002602	TITANIUM PLATE, 5- HOLES, BROAD MEDIUM, PALMAR, LEFT, SYSTEM 2.5	1	15
48	4002603	TITANIUM PLATE, 5 - HOLES, BROAD MEDIUM, PALMAR, RIGHT, SYSTEM 2.5	1	15
49	4002604	TITANIUM PLATE, 6 - HOLES, BROAD LARGE, PALMAR, LEFT, SYSTEM 2.5	1	15
50	4002605	TITANIUM PLATE, 6 - HOLES, BROAD LARGE, PALMAR, RIGHT, SYSTEM 2.5	1	15
51	4002606	TITANIUM LONG XL PLATE, 11 - HOLES NARROW, PALMAR, LEFT, SYSTEM 2.5	1	17
52	4002607	TITANIUM LONG XL PLATE, 11 - HOLES NARROW, PALMAR, RIGHT, SYSTEM 2.5	1	17
53	4002608	TITANIUM PLATE, 2/4-HOLES, L-FORM, LEFT, SYSTEM 2.5	1	18
54	4002609	TITANIUM PLATE, 2/4-HOLES, L-FORM, RIGHT, SYSTEM 2.5	1	18
55	4002610	TITANIUM PLATE, 6-HOLES, STRAIGHT, SYSTEM 2.5	1	18
56	4002611	TITANIUM PLATE, 8-HOLES, T-FORM, SYSTEM 2.5	1	19
57	4002612	TITANIUM PLATE, 6-HOLES, HOOK, SYSTEM 2.5	1	19
58	4002613	TITANIUM LONG L PLATE, 7 - HOLES NARROW, PALMAR, RIGHT, SYSTEM 2.5	1	16
59	4002614	TITANIUM LONG L PLATE, 7 - HOLES NARROW, PALMAR, LEFT, SYSTEM 2.5	1	16
60	4002615	TITANIUM ATRAUMATIC FRAME PLATE, 14-HOLES, SIZE 1, RIGHT, SYSTEM 2.5	1	20
61	4002616	TITANIUM ATRAUMATIC FRAME PLATE, 14-HOLES, SIZE 1, LEFT, SYSTEM 2.5	1	20
62	4002617	TITANIUM ATRAUMATIC FRAME PLATE, 14-HOLES, SIZE 2, RIGHT, SYSTEM 2.5	1	20
63	4002618	TITANIUM ATRAUMATIC FRAME PLATE, 14-HOLES, SIZE 2, LEFT, SYSTEM 2.5	1	20
Screws				
64	5028010	TITANIUM LOCKING SCREW, 2,5 X 10 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
65	5028012	TITANIUM LOCKING SCREW, 2,5 X 12 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
66	5028014	TITANIUM LOCKING SCREW, 2,5 X 14 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
67	5028016	TITANIUM LOCKING SCREW, 2,5 X 16 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21

Pos.	Cat.No.	Description	Qty.	Picture page
68	5028018	TITANIUM LOCKING SCREW, 2,5 X 18 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
69	5028020	TITANIUM LOCKING SCREW, 2,5 X 20 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
70	5028022	TITANIUM LOCKING SCREW, 2,5 X 22 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
71	5028024	TITANIUM LOCKING SCREW, 2,5 X 24 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
72	5028026	TITANIUM LOCKING SCREW, 2,5 X 26 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
73	5028028	TITANIUM LOCKING SCREW, 2,5 X 28 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
74	5028030	TITANIUM LOCKING SCREW, 2,5 X 30 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
75	5028110	TITANIUM SCREW STANDARD 2,5 X 10 MM SELF-TAPPING, TX 7, PACK/5	1	21
76	5028112	TITANIUM SCREW STANDARD 2,5 X 12 MM SELF-TAPPING, TX 7, PACK/5	2	21
77	5028114	TITANIUM SCREW STANDARD 2,5 X 14 MM SELF-TAPPING, TX 7, PACK/5	2	21
78	5028116	TITANIUM SCREW STANDARD 2,5 X 16 MM SELF-TAPPING, TX 7, PACK/5	2	21
79	5028118	TITANIUM SCREW STANDARD 2,5 X 18 MM SELF-TAPPING, TX 7, PACK/5	2	21
80	5028120	TITANIUM SCREW STANDARD 2,5 X 20 MM SELF-TAPPING, TX 7, PACK/5	1	21
81	5028122	TITANIUM SCREW STANDARD 2,5 X 22 MM SELF-TAPPING, TX 7, PACK/5	1	21
82	5028124	TITANIUM SCREW STANDARD 2,5 X 24 MM SELF-TAPPING, TX 7, PACK/5	1	21
83	5028126	TITANIUM SCREW STANDARD 2,5 X 26 MM SELF-TAPPING, TX 7, PACK/5	1	21
84	5028128	TITANIUM SCREW STANDARD 2,5 X 28 MM SELF-TAPPING, TX 7, PACK/5	1	21
85	5028130	TITANIUM SCREW STANDARD 2,5 X 30 MM SELF-TAPPING, TX 7, PACK/5	1	21
86	5028210	TITANIUM SCREW 2,5 X 10 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
87	5028212	TITANIUM SCREW 2,5 X 12 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
88	5028214	TITANIUM SCREW 2,5 X 14 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
89	5028216	TITANIUM SCREW 2,5 X 16 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
90	5028218	TITANIUM SCREW 2,5 X 18 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
91	5028220	TITANIUM SCREW 2,5 X 20 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
92	5028222	TITANIUM SCREW 2,5 X 22 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
93	5028224	TITANIUM SCREW 2,5 X 24 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
94	5028226	TITANIUM SCREW 2,5 X 26 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
95	5028228	TITANIUM SCREW 2,5 X 28 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
96	5028230	TITANIUM SCREW 2,5 X 30 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22

S1486 - SMALL KING KONG SYSTEM 2.5 SET



Pos.	Cat.No.	Description	Qty.	Picture page
Containers				
1	1301037	AGOPAQX2 TRAY FOR IMPLANT MODULES W/ 2 LIDS, W/O IMPLANT MODULES	1	45
2	1000379	AGOPAQX4 PLATE-TRAY SYSTEM 2.5, NARROW, W/O IMPLANTS	1	47
3	1000438	AGOPAQX4 SCREW TRAY SYSTEM 2.5, W/O IMPLANTS	1	46
4	1000381	AGOPAQX4 PLATE-TRAY SYSTEM 2.5, W/O IMPLANTS	1	46
5	1000384	AGOPAQX4 SCREW TRAY SYSTEM 2.5, W/O IMPLANTS	1	46
Instruments				
6	1001910	SCREWDRIVER BLADE, DOUBLE TX 7 NON-CANNULATED, 90 MM, AO-SHAFT	1	48
7	1001911	SCREWDRIVER BLADE, DOUBLE TX 7 CANNULATED, 90 MM, AO-SHAFT	1	48
8	1004036	SILICONE SCREWDRIVERHANDLE BLACK CANNULATED, 12 CM, AO-SHAFT	1	49
9	1002527	PLATE- /SCREWHOLDING FORCEPS ANGLED, 15CM/6"	1	49
10	1006024	AGOMED DEPTH GAUGE 40 MM, ROUND F. 2,0 - 4,5 MM SCREWS	1	49
11	1008009	HANDLE FOR DRILL-GUIDE, ONLY HAND-SYSTEMS	1	49
12	1008040	DRILL-GUIDE ONLY, VARIABLE, LENGTH 35 MM, DRILL Ø 2.0, ROUND TIP	1	51
13	1008050	K-WIRE GUIDE ONLY, VARIABLE, K-WIRE Ø 1.1, LENGTH 35 MM, ROUND TIP	1	51
14	1008051	K-WIRE GUIDE ONLY, SCREWABLE, K-WIRE Ø 1.1, LENGTH 45 MM	1	51
15	1008054	AIMING DEVICE RADIUS 2.5, NARROW, LEFT	1	52
16	1008055	AIMING DEVICE RADIUS 2.5, NARROW, RIGHT	1	52

Pos.	Cat.No.	Description	Qty.	Picture page
17	1008065	DRILL-GUIDE ONLY, SCREWABLE, LENGTH 45 MM, DRILL Ø 2.0	1	51
18	1008041	DRILL GUIDE ONLY, SCREWABLE, DRILL Ø2.0, 0° TO 10°, SYSTEM 2.5	1	50
19	1008066	DRILL-GUIDE ONLY, VARIABLE, LENGTH 35 MM, DRILL Ø 2.0	1	51
20	1008068	K-WIRE GUIDE ONLY, VARIABLE, K-WIRE Ø 1.1, LENGTH 35 MM	1	51
21	1008069	SLEEVE REDUCTION GUIDE, 35 MM LENGTH, K-WIRE Ø1,1, SYSTEM 2.5	1	50
22	1101340	OLIVE KIRSCHNER-WIRES, THREAD 10 MM, TROCAR/ROUND, 60 X 1.1 MM, SYSTEM 2.5, PACK/1	3	50
23	1019106	TWIST DRILL 2,0 X 95 MM, CANNULATED 1,2 MM, AO-SHAFT, STOP 26 MM	1	48
24	1101023	KIRSCHNER WIRE 100 X 1,1 MM TROCAR/ROUND, PACK/6	1	50
25	1201008	TWIST DRILL 2,0 X 103 MM, AO-SHAFT, GREY CODE	1	48
26	1101063	WERBER CBS-K-WIRE 150 X 1,4 MM TROCAR / ROUND, PACK/6	1	50
27	1104009	AGOMED K-WIRE DISPENSER 15 CM/6" FOR K-WIRE 0,8 - 1,6 MM DIAMETER	2	50
Plates				
28	4002594	TITANIUM PLATE, 4-HOLES, NARROW SMALL, PALMAR, LEFT, SYSTEM 2.5	1	12
29	4002595	TITANIUM PLATE, 4-HOLES, NARROW SMALL, PALMAR, RIGHT, SYSTEM 2.5	1	12
30	4002596	TITANIUM PLATE, 5-HOLES, NARROW MEDIUM, PALMAR, LEFT, SYSTEM 2.5	1	13
31	4002597	TITANIUM PLATE, 5-HOLES, NARROW MEDIUM, PALMAR, RIGHT, SYSTEM 2.5	1	13
32	4002598	TITANIUM PLATE, 6-HOLES, NARROW LARGE, PALMAR, LEFT, SYSTEM 2.5	1	13
33	4002599	TITANIUM PLATE, 6-HOLES, NARROW LARGE, PALMAR, RIGHT, SYSTEM 2.5	1	13



Pos.	Cat.No.	Description	Qty.	Picture page
34	4002608	TITANIUM PLATE, 2/4-HOLES, L-FORM, LEFT, SYSTEM 2.5	1	18
35	4002609	TITANIUM PLATE, 2/4-HOLES, L-FORM, RIGHT, SYSTEM 2.5	1	18
36	4002610	TITANIUM PLATE, 6-HOLES, STRAIGHT, SYSTEM 2.5	1	18
37	4002611	TITANIUM PLATE, 8-HOLES, T-FORM, SYSTEM 2.5	1	19
38	4002612	TITANIUM PLATE, 6-HOLES, HOOK, SYSTEM 2.5	1	19
Screws				
39	5028010	TITANIUM LOCKING SCREW, 2,5 X 10 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
40	5028012	TITANIUM LOCKING SCREW, 2,5 X 12 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
41	5028014	TITANIUM LOCKING SCREW, 2,5 X 14 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
42	5028016	TITANIUM LOCKING SCREW, 2,5 X 16 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
43	5028018	TITANIUM LOCKING SCREW, 2,5 X 18 MM MULTIDIRECTIONAL, TX 7, PACK/5	2	21
44	5028020	TITANIUM LOCKING SCREW, 2,5 X 20 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
45	5028022	TITANIUM LOCKING SCREW, 2,5 X 22 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
46	5028024	TITANIUM LOCKING SCREW, 2,5 X 24 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
47	5028026	TITANIUM LOCKING SCREW, 2,5 X 26 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
48	5028028	TITANIUM LOCKING SCREW, 2,5 X 28 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
49	5028030	TITANIUM LOCKING SCREW, 2,5 X 30 MM MULTIDIRECTIONAL, TX 7, PACK/5	1	21
50	5028110	TITANIUM SCREW STANDARD 2,5 X 10 MM SELF-TAPPING, TX 7, PACK/5	1	21
51	5028112	TITANIUM SCREW STANDARD 2,5 X 12 MM SELF-TAPPING, TX 7, PACK/5	2	21
52	5028114	TITANIUM SCREW STANDARD 2,5 X 14 MM SELF-TAPPING, TX 7, PACK/5	2	21
53	5028116	TITANIUM SCREW STANDARD 2,5 X 16 MM SELF-TAPPING, TX 7, PACK/5	2	21
54	5028118	TITANIUM SCREW STANDARD 2,5 X 18 MM SELF-TAPPING, TX 7, PACK/5	2	21

Pos.	Cat.No.	Description	Qty.	Picture page
55	5028120	TITANIUM SCREW STANDARD 2,5 X 20 MM SELF-TAPPING, TX 7, PACK/5	1	21
56	5028122	TITANIUM SCREW STANDARD 2,5 X 22 MM SELF-TAPPING, TX 7, PACK/5	1	21
57	5028124	TITANIUM SCREW STANDARD 2,5 X 24 MM SELF-TAPPING, TX 7, PACK/5	1	21
58	5028126	TITANIUM SCREW STANDARD 2,5 X 26 MM SELF-TAPPING, TX 7, PACK/5	1	21
59	5028128	TITANIUM SCREW STANDARD 2,5 X 28 MM SELF-TAPPING, TX 7, PACK/5	1	21
60	5028130	TITANIUM SCREW STANDARD 2,5 X 30 MM SELF-TAPPING, TX 7, PACK/5	1	21
61	5028210	TITANIUM SCREW 2,5 X 10 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
62	5028212	TITANIUM SCREW 2,5 X 12 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
63	5028214	TITANIUM SCREW 2,5 X 14 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
64	5028216	TITANIUM SCREW 2,5 X 16 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
65	5028218	TITANIUM SCREW 2,5 X 18 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	2	22
66	5028220	TITANIUM SCREW 2,5 X 20 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
67	5028222	TITANIUM SCREW 2,5 X 22 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
68	5028224	TITANIUM SCREW 2,5 X 24 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
69	5028226	TITANIUM SCREW 2,5 X 26 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
70	5028228	TITANIUM SCREW 2,5 X 28 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22
71	5028230	TITANIUM SCREW 2,5 X 30 MM CANNULATED, LOCKING, FULLTHREAD, TX 7, PACK/5	1	22

MORE OF OUR PRODUCTS

2010015

Cannulated Screws




2010033

AgoFix 4.0 system



2010003

Titanium plating system for craniomaxillofacial osteosynthesis



2010002

Foot & Ankle system



2010030

Pediatric Screw system 6.2



2010008

AgoFix Neuro Plating system 1.5




NEW 2010053

Anterior 4.0 system



2010017

AGO-Pediatric system 3.5 / 4.0



2010016

Titanium Mesh systems for craniomaxillofacial osteosynthesis / Neuro





AGOMED MEDIZIN-TECHNIK GMBH
ÖSCHWEG 29 - D-78567 FRIDINGEN - GERMANY

☎ **(+49) 7463 / 267 06 16**
✉ **info@agomed.com**
🌐 **<http://www.agomed.com>**

MADE IN GERMANY

Ref.num. 2010065
Revision 01 / 12.2020

