



FOOT & ANKLE STRAIGHT PLATES

EN Working with AGOMED means working with a company committed to excellence. Our products are German engineered, biocompatible, and cutting-edge. Our team members collectively bring decades of medical technology experience to the table.

We're reliable. We're flexible. We're inventive.

AGOMED collaborates with distributors, hospitals, and doctors worldwide to create and develop exacting, state-of-the-art solutions. We work closely with specialists to insure that our implants for trauma and extremities improve the patient's quality of life.

Patient safety is always our number one priority.

DE Mit AGOMED zu arbeiten, bedeutet mit einem Unternehmen zu arbeiten, dass sich zur Exzellenz und Hochleistung verpflichtet.

Unsere Produkte werden von deutschen Ingenieuren entwickelt, sind biokompatibel, innovativ und auf dem neuesten Stand der Technik. Unser Team hat jahrzehntelange Erfahrung im Bereich der Medizintechnik.

Wir sind zuverlässig. Wir sind flexibel. Wir sind innovativ.

AGOMED arbeitet weltweit mit Unternehmen, Krankenhäusern und Ärzten zusammen, um anspruchsvolle, hochmoderne Lösungen zu entwickeln. Wir arbeiten eng mit Spezialisten zusammen, um sicherzustellen, dass unsere Implantate für Traumatologie und Extremitäten die Lebensqualität der Patienten verbessern.

Die Sicherheit der Patienten steht bei uns immer an erster Stelle.

ES Trabajar con AGOMED significa trabajar con una empresa comprometida con la excelencia. Nuestros productos son diseñados en Alemania, biocompatibles y de vanguardia. Los miembros de nuestro equipo aportan, de manera colectiva, décadas de experiencia en tecnología médica.

Somos fiables. Somos flexibles. Somos innovadores.

AGOMED colabora con distribuidores, hospitales y doctores a nivel mundial con el objetivo de crear y desarrollar soluciones rigurosas y con lo último de la tecnología. Cooperamos estrechamente con especialistas para garantizar que nuestros implantes para trauma y extremidades mejoren la calidad de vida de los pacientes.

La seguridad de los pacientes siempre es nuestra prioridad.

QUALITY IS PART OF OUR DNA



DIN EN ISO 13485

AGOMED MEANS QUALITY. WE EXCEED ISO 13485 STANDARDS AND HAVE CE CERTIFICATION.

AGOMED BEDEUTET QUALITÄT. WIR SIND ZERTIFIZIERT NACH ISO 13485 UND HABEN DIE CE ZERTIFIZIERUNG

AGOMED ES SINÓNIMO DE CALIDAD. SUPERAMOS LOS ESTÁNDARES ISO 13485 Y CONTAMOS CON CERTIFICACIÓN CE.



CE CERTIFICATE

INDEX / INDEX / ÍNDICE

STRAIGHT FOOT PLATING SYSTEM 2.7

GERADES FUSS PLATTEN SYSTEM 2.7
SISTEMA RECTO DE PLACAS DE PIE 2.7

2

TITANIUM SCREWS 2.7

TITANSCHRAUBEN 2.7
TORNILLOS DE TITANIO 2.7

8

INDEX OF ITEM NUMBERS

INDEX VON ARTIKEL NUMMERN
ÍNDICE DE NÚMEROS DE ARTICULO

10



EVANS TITANIUM PLATE /
EVANS TITAN PLATTE / PLACA DE TITANIO EVANS



SAMPLE INDICATIONS:

- Osteotomy
- Arthrodesis
- Joint fusion
- Fracture repair
- Fracture fixation
- Ligament fixation
- Bone reconstruction

1.1 mm



4004378
2-holes, small size
2-loch, groesse S
2 agujeros, pequeña



4004379
2-holes, medium size
2-loch, groesse M
2 agujeros, intermedio



4004380
2-holes, large size
2-loch, groesse L
2 agujeros, grande

EVANS TITANIUM COMPRESSION PLATE /
EVANS KOMPRESSIONS TITAN PLATTE / PLACA DE TITANIO EVANS CON COMPRESIÓN

1.1 mm

COMPRESSION HOLES




4004381
2-holes, compression, small size
2-loch, kompressions, groesse S
2 agujeros, compresión, pequeña



4004384
3-holes, compression, small size
3-loch, kompressions, groesse S
3 agujeros, compresión, pequeña



4004382
2-holes, compression, medium size
2-loch, kompressions, groesse M
2 agujeros, compresión, intermedio



4004385
3-holes, compression, medium size
3-loch, kompressions, groesse M
3 agujeros, compresión, intermedio



4004383
2-holes, compression, large size
2-loch, kompressions, groesse L
2 agujeros, compresión, grande




4004386
3-holes, compression, large size
3-loch, kompressions, groesse L
3 agujeros, compresión, grande



AGOSALIS TITANIUM PLATE /
AGOSALIS TITAN PLATTE / PLACA DE TITANIO AGOSALIS

SAMPLE INDICATIONS:

- Interposition and revision arthrodesis of Lisfranc joints
- Iliac crest interposition arthrodesis of first toe
- Axial corrections and metatarsal shortening
- Revision arthrodesis of the MTP Joint
- Arthrodesis of first MTP Joint
- ORIF in trauma cases
- Lisfranc arthrodesis



STRAIGHT / GERADE / RECTA $\overline{1.6\text{ mm}}$



4001086
4-holes
4-loch
4-agujeros




4001087
5-holes
5-loch
5-agujeros




AGOSALIS TITANIUM PLATE /
AGOSALIS TITAN PLATTE / PLACA DE TITANIO AGOSALIS

SAMPLE INDICATIONS:


- Multiple Arthrodesis from MT to Os Cuneiforme



STRAIGHT / GERADE / RECTA $\overline{1.6\text{ mm}}$



4001101
10-holes, pack/1
10-loch, pack/1
10-agujeros, paq./1



4001102
12-holes, pack/1
12-loch, pack/1
12-agujeros, paq./1



4001103
14-holes, pack/1
14-loch, pack/1
14-agujeros, paq./1



2.7 TITANIUM PLATE /
2.7 TITAN PLATTE / PLACA DE TITANIO 2.7



SAMPLE INDICATION:
• Lisfranc Joint

L FORM $\frac{1.6 \text{ mm}}$




STRAIGHT / GERADE / RECTA $\frac{1.6 \text{ mm}}$

T FORM $\frac{1.6 \text{ mm}}$



STRAIGHT / GERADE / RECTA $\frac{2.0 \text{ mm}}$






TITANIUM SCREWS 2.7 /
TITANSCHRAUBEN 2.7 / TORNILLOS DE TITANIO 2.7



L	Ø 2.7 mm	Ø 2.7 mm	Pack/Pack/Paq.
8 mm	5001108	5001008	5
10 mm	5001110	5001010	5
12 mm	5001112	5001012	5
14 mm	5001114	5001014	5
16 mm	5001116	5001016	5
18 mm	5001118	5001018	5
20 mm	5001120	5001020	5
22 mm	5001122	5001022	5
24 mm	5001124	5001024	5
26 mm	5001126	5001026	5
28 mm	5001128	5001028	5
30 mm	5001130	5001030	5




MORE OF OUR PRODUCTS /
MEHR ÜBER UNSERE PRODUKTE / MÁS DE NUESTROS PRODUCTOS

2010002
Foot & Ankle
system



2010015
Cannulated
Screws



2010003
Titanium plating system
for craniomaxillofacial osteosynthesis



2010030
Pediatric
Screw system 6.2




2010004
Plating system for
Upper Extremities



2010008
AgoFix Neuro
Plating system 1.5



2010017
AGO-Pediatric
system 3.5 / 4.0



2010007
King Kong Multidirectional
Radius system 2.5



2010016
Titanium Mesh systems for
craniomaxillofacial osteosynthesis / Neuro





NOTES / NOTIZEN / NOTAS

Cat.No.	Page
5001110	8
5001112	8
5001114	8
5001116	8
5001118	8
5001120	8
5001122	8
5001124	8
5001126	8
5001128	8
5001130	8

NOTES / NOTIZEN / NOTAS

NOTES / NOTIZEN / NOTAS







NOTES / NOTIZEN / NOTAS

Lined area for notes on page 12, featuring horizontal blue lines and a background pattern of faint, stylized 'A' logos.



NOTES / NOTIZEN / NOTAS

Lined area for notes on page 13, featuring horizontal blue lines and a background pattern of faint, stylized 'A' logos.



AGOMED MEDIZIN-TECHNIK GMBH
ÖSCHWEG 29 - D-78567 FRIDINGEN - GERMANY

☎ (+49) 7463 / 267 06 16

✉ info@agomed.com

🌐 <http://www.agomed.com>

MADE IN GERMANY

Ref.num. 2010047
Revision 01 / 05.2020

